

■■ NavTools: Assistive Gaze Tracking & Voice Assistant

A Hands-Free Assistive Control Suite using MediaPipe Gaze Tracking and Offline Voice Command Orchestration

Group No. 7 · 8th Semester Major Project · 2026

■ Table of Contents

- Overview
 - System Architecture
 - Core Modules
 - Gaze Tracking and Clicking Gestures
 - Jim Interactive Voice Assistant
 - Animated Orb UI
 - General Settings Tab
 - Software Setup
 - How to Run
 - Quick Start (run.bat)
 - Multimodal Console Launcher
 - Control Center Dashboard
 - Unified State Management
 - Voice Assistant Reference Guide
 - Troubleshooting
 - License
-

■ Overview

NavTools is a completely non-invasive, premium assistive control suite that empowers individuals with physical impairments to operate computers hands-free. The system merges **webcam-based eye gaze tracking** with an **offline voice-activated assistant (Jim)**.

Unlike legacy biosignal-based interfaces (like EEG/EOG), NavTools requires **zero extra hardware**. It leverages standard computer webcams and microphones to deliver smooth cursor movement, rich facial-gesture-based clicking, and highly responsive voice command execution.

Key Highlights



■ Core Modules

1. Camera-Based Gaze Tracking

The **Gaze Tracker** acts as a virtual mouse by tracking the user's face and eyes using MediaPipe's High-Fidelity Landmark Model (`face_landmarker.task`).

- **Exponential Moving Average (EMA) Smoothing:** Eliminates natural micro-saccades and camera jitter to provide smooth, organic cursor movement.
- **Eye Aspect Ratio (EAR) Clicking:** Detects natural blinks and winks for full mouse emulation:
 - ■ **Left Eye Wink:** Executes a **Left Mouse Click**.
 - ■ **Right Eye Wink:** Executes a **Right Mouse Click**.
 - ■ ■ **Double Blink:** Executes a **Double-Click**.
 - ■ ■ **Extended Hold Blink:** Triggers mouse **Click & Drag**.

2. Jim Interactive Voice Assistant

Jim is an offline-first interactive digital companion that processes vocal commands. It features:

- **120+ Voice Commands:** Spanning app launches, browser/tab navigation, file operations, system controls, YouTube Music, window management, and more.
- **Smart Calculation Parser:** Extracts dynamic spoken expressions (e.g., *"calculate fifteen times three plus four"*) and executes them safely.
- **Smart Web & Local Search:** Supports targeted phrases (e.g., *"search Google for python guides"*, *"find file design_specifications"*, *"type hello world"*).
- **Positional Link/File Clicking:** Commands like *"click first result"*, *"click second link"*, *"open third file"* work in browsers and File Explorer.
- **Kokoro TTS Engine:** Integrates lightweight ONNX models to produce crystal-clear human-like speech offline.
- **Dynamic Mic Sensitivity:** Reads `attention.mic_sensitivity` before each listen call for real-time adjustment.

3. Animated Orb UI (PyQt5 + QWebEngineView)

The system renders a **floating, frameless, fully transparent orb** on your screen using PyQt5's `QMainWindow` + `QWebEngineView`:

- **Frameless & Transparent:** `Qt.FramelessWindowHint` | `Qt.WA_TranslucentBackground` — no border, pure floating orb.
- **Always On Top:** `Qt.WindowStaysOnTopHint` — remains visible above all other windows.
- **Drag-to-Move:** Click and drag the orb to any screen position; delta coordinates are relayed via `console.log("ORB_DRAG:dx,dy")` to the Python bridge.

- **Color-coded States:**
 - ■ **Idle** — Dim pulsing blue (waiting for wake word)
 - ■ **Listening** — Bright cyan/green pulse (hearing your command)
 - ■ **Speaking** — Amber glow (Jim is responding)
- **Double-click** → Opens the Tkinter Settings Dashboard.
- **Right-click** → Exits the entire application gracefully.

4. General Settings Tab

A new **General Settings** tab in the Control Center Dashboard provides runtime-adjustable parameters:

Setting	Description	Default
Voice Profile	Selects TTS voice: "Soft Female (Zira)" or "Natural Female (Sarah)"	Soft Female (Zira)
Speaking Speed	Adjusts voice speed from 0.5× to 2.0×	1.0×
Microphone Sensitivity	Energy threshold for voice pick-up (50–1000); lower = more sensitive	300

Changes apply instantly — the voice assistant reads these from **AttentionState** before each speak/listen cycle.

■ Software Setup

Prerequisites

- **OS:** Windows 10/11 (preferred for pywin32 API system hooks and SAPI5 TTS).
- **Python:** 3.9 through 3.11.
- **Webcam & Microphone:** Any standard built-in or external USB device.

Installation

Option A — Automated Installer (Recommended)

Simply double-click **install.bat**. It will:

- Check for Python in PATH.
- Create a **.venv** virtual environment.
- Upgrade pip and install all dependencies from **requirements.txt**.
- Download the MediaPipe Face Landmarker model (~80MB) into **data/**.

Option B — Manual Setup

- **Clone the Repository:**

```
```bash
```

```
git clone https://github.com/JatinderpalSingh9321/Assistant-Project-for-Paralysed-people.git
```

```
cd Assistant-Project-for-Paralysed-people
```

```
```
```

- **Set Up Virtual Environment:**

```
```bash
```

```
python -m venv .venv
```

```
.\.venv\Scripts\activate
```

```
```
```

- **Install Dependencies:**

```
```bash
```

```
pip install -r requirements.txt
```

```
```
```

- **Download Face Landmarker Model:**

```
```bash
```

```
PowerShell
```

```
Invoke-WebRequest -Uri "https://storage.googleapis.com/mediapipe-models/face_landmarker/face_landmarker/float16/1/face_landmarker.task" -OutFile "data\face_landmarker.task"
```

```
```
```

- *(Optional)* **Download Kokoro TTS Models:**

```
```bash
```

```
python scripts/download_kokoro.py
```

```
```
```

■ How to Run

Method 1: Quick Start (run.bat)

Double-click `run.bat` — it automatically uses the `.venv` Python and launches the full Control Center dashboard:

```
run.bat
```

Method 2: Standalone Multimodal Console Launcher

Runs both Eye Gaze Tracker and Voice Assistant as independent threads with console logs:

```
# Basic start (Gaze + Voice)
python -m src.multimodal_launcher

# Start with a real-time webcam preview panel for camera calibration
python -m src.multimodal_launcher --preview

# Disable gaze module (voice only)
python -m src.multimodal_launcher --no-gaze
```

Method 3: Control Center Dashboard

Launches the full dark-themed GUI control panel with three settings tabs:

```
python -m src.gui_app
```

The dashboard includes:

- **Voice Assistant tab** — Start/Stop Jim, enable Attention Gating.
- **Gaze Tracker tab** — Configure camera index, EMA smoothing, toggle live preview.
- **General Settings tab** — Select voice profile, speaking speed, microphone sensitivity.

■ Unified State Management

To prevent cursor drift or voice triggers when the user is not actively looking at the screen, NavTools implements a thread-safe singleton state class, **AttentionState**.

AttentionState Properties

| Property | Type | Description |
|------------------------------|--------------------|--|
| <code>is_attentive</code> | <code>bool</code> | Whether the user is focused on the screen (staleness check after 2s) |
| <code>gaze_x / gaze_y</code> | <code>float</code> | Normalized gaze coordinates (0.0–1.0) |
| <code>voice_name</code> | <code>str</code> | Active TTS voice profile (" zira " or " sarah ") |
| <code>voice_speed</code> | <code>float</code> | Speaking speed multiplier (0.5–2.0) |
| <code>mic_sensitivity</code> | <code>int</code> | Recognizer energy threshold (50–1000) |

If **Attention Gating** is enabled:

- **Attentive Calibration:** The eye gaze tracker checks if the user's pupils are focused on the monitor.

- **Signal Gating:** If the user looks away, `attention.is_attentive` is set to `False`.
- **Trigger Lockout:** The Voice Assistant pauses microphone streaming and cursor positioning freezes.
- **Resumption:** Simply looking back at the screen immediately resumes all controls.

■■ Voice Assistant Reference Guide

Jim listens for the wake phrase **"wake up Jim"** or **"hey Jim"**. Once activated, speak any of the following command classes:

Browser & Tab Navigation

| Command | Action |
|---|------------------------------------|
| <i>"open google" / "open youtube"</i> | Launches browser to specified URLs |
| <i>"new tab" / "close tab"</i> | Manages browser tabs |
| <i>"next tab" / "previous tab"</i> | Cycles through active tabs |
| <i>"go back" / "go forward" / "refresh page"</i> | Standard page navigation |
| <i>"click first result" / "click second link"</i> | Clicks nth link on current page |

Local Application Launchers

- *"open calculator" / "open notepad" / "open task manager" / "open terminal"*
- *"open windows settings" / "open device manager" / "open remote desktop"*
- *"open word" / "open excel" / "open powerpoint" / "open outlook"*
- *"open steam" / "open discord" / "open spotify" / "open vs code" / "open photoshop"*
- *"open downloads" / "open documents" / "open pictures" / "open desktop"*

Windows & System Controls

- *"minimize window" / "maximize window" / "snap left" / "snap right" / "full screen"*
- *"scroll down" / "scroll up" / "page up" / "page down" / "go to top" / "go to bottom"*
- *"take screenshot"* — Windows Snip & Sketch overlay.
- *"lock screen"* — Instantly locks the Windows session.
- *"volume up" / "volume down" / "mute"* — System-wide audio control.
- *"copy" / "paste" / "undo" / "redo" / "select all"*

Global Media Controls

- *"play" / "pause" / "play pause"* — Windows global media keys.
- *"next song" / "previous song" / "skip track"* — Track navigation.

YouTube Music Controls

| Command | Action |
|---|--------------------------------------|
| "play music" / "pause music" | Toggle playback on YouTube Music tab |
| "next music" / "previous music" | Skip tracks |
| "mute music" / "unmute music" | Mute/unmute |
| "like this song" / "dislike this song" | Rate current track |
| "shuffle music" / "shuffle on" | Toggle shuffle |
| "volume up music" / "volume down music" | Adjust music volume |
| "open youtube music" | Opens music.youtube.com in browser |

File Explorer Controls

- "open first file" / "open second file" ... "open fifth file" / "open last file"
- "select first file" / "next file" / "previous file"
- "rename file" / "delete file" / "copy file" / "paste file" / "new folder"

Smart Dynamic Handlers

- "calculate [expression]" — Solves equations vocally.
- "search for [query]" — Smart search (Google or YouTube Music depending on context).
- "google search for [query]" — Always searches Google.
- "search in file explorer for [name]" — Opens File Explorer with query.
- "find file [filename]" — File Explorer structured search.
- "type [text]" — Types text directly into the active cursor field.
- "open [app name]" / "close [app name]" — Dynamic app launch/close.
- "play [song/artist name]" — Searches YouTube Music.

Gaze Tracker Controls (via Voice)

- "launch eye tracking" / "start eye cursor" — Starts the gaze tracking module.
- "stop eye tracking" / "close eye cursor" — Stops the gaze tracking module.

Settings & Assistant Controls

- "open settings" / "show dashboard" / "open control center" — Opens the Control Center.
- "close settings" / "hide panel" — Hides the Control Center.
- "help" / "what can you do" — Jim lists available commands.
- "stop listening" / "go to sleep" — Puts Jim to sleep (wake phrase restores).
- "close the assistant" / "exit application" — Graceful full shutdown.

■ Troubleshooting

| Issue | Likely Cause | Solution |
|------------------------|--|---|
| Jittery cursor | Low room lighting or low camera FPS | Increase ambient lighting; increase gaze SMOOTHING value (e.g., 0.9) in the Gaze Tracker settings tab. |
| Blinks not clicking | Face landmark model file is missing | Ensure data/face_landmarker.task is present. Run install.bat or download manually. |
| Voice command lag | High ambient microphone noise | Move to a quieter area or raise Microphone Sensitivity in General Settings. |
| Jim not waking up | Energy threshold too high | Lower Microphone Sensitivity in the General Settings tab (try 150–250). |
| Kokoro TTS unavailable | Model weights missing or not installed | System auto-falls back to Windows SAPI5. Run python scripts/download_kokoro.py to install Kokoro. |
| Orb not appearing | PyQt5/QWebEngine not installed | Run pip install PyQt5 PyQtWebEngine or re-run install.bat . |
| run.bat fails | .venv not created | Run install.bat first, or create .venv manually and install requirements.txt . |

■ Project Structure

Assistant-Project-for-Paralysed-people/

```

■■■■ src/
■   ■■■■ gui_app.py           # Main app: PyQt5 Orb + Tkinter Dashboard
■   ■■■■ voice_assistant.py  # Jim: 120+ voice commands, Kokoro/SAPI5 TTS
■   ■■■■ attention_state.py  # Thread-safe shared state singleton
■   ■■■■ gaze_tracker.py     # MediaPipe face landmark gaze + EAR clicking
■   ■■■■ multimodal_launcher.py # CLI launcher (gaze + voice threads)
■   ■■■■ eye_tracker.py      # Low-level eye tracking utilities
■   ■■■■ gaze_mouse.py       # PyAutoGUI mouse control bridge
■   ■■■■ utils.py            # Logger and path utilities
■■■■ scripts/
■   ■■■■ generate_report.py   # B.Tech major project report generator (DOCX)
■   ■■■■ download_kokoro.py   # Kokoro ONNX model downloader
■   ■■■■ test_kokoro.py       # Kokoro TTS test script
■■■■ data/
■   ■■■■ face_landmarker.task # MediaPipe model (downloaded by install.bat)
■■■■ models/
■   ■■■■ kokoro/              # Kokoro ONNX model weights (optional)
■■■■ orb.html                 # Animated orb HTML/CSS/JS
■■■■ requirements.txt          # Python dependencies
■■■■ install.bat              # One-click installer

```

```
■■■ run.bat # One-click launcher (uses .venv)
■■■ BCI_NavTools_Major_Project_Report.docx # Full B.Tech project report
```

■ License

This project is licensed under the MIT License.

Assistive Technology made elegant, lightweight, and accessible to everyone. ■■■■

Group 7 — 8th Semester Major Project, 2026